

**Bring your sample in the same day it is drawn.  
Return to Dixie County Health Department on  
Monday or Thursday before 12:00 p.m. (Noon).  
Fee is \$20.00 for each sample**

Name: \_\_\_\_\_ Phone: \_\_\_\_\_

Mailing Address: \_\_\_\_\_

Property Location/Subdivision: \_\_\_\_\_

Date and Time of Collection: \_\_\_\_\_

Type of Supply: Private \_\_\_\_\_ Well Bottle \_\_\_\_\_ Other \_\_\_\_\_

Type of Sample: Compliance \_\_\_\_\_ Well Survey \_\_\_\_\_

Check One: \_\_\_\_\_ Raw \_\_\_\_\_ Distribution


When collecting a water sample:

1. Sample your well water from a point used frequently.
2. Run water for 10 to 15 minutes before collecting.
3. Follow the instructions for sample bag below. **Be sure to fill water to top fill line.**


**Instructions for Use of Nasco's Sodium Thiosulfate Whirl-Pak<sup>®</sup> Bags**

*Caution: Do not allow the thio tablet to fall out of the bag.*

**OPENING**




1. Tear off top at scored line.




2. Pull tabs outward to open bag. Sometimes a pull on the bottom is also helpful.

**FILLING**




3. Fill bag to the fill line.


**CLOSING**



4. Pull wire ends to close bag.



5. Whirl bag 3 complete revolutions. (Do not roll tapes down.) Leave airspace for mixing in bag.



6. Turn top wire inward on opposite face of fold.

**NASCO, FORT ATKINSON, WISCONSIN 53538 • MODESTO, CALIFORNIA 95352**  
U 5352 PRINTED IN U.S.A.

Dixie County Health Department  
Phone (352) 498-1360 • 1530 SE 12<sup>th</sup> Avenue • Cross City, FL 32628